

Allergies and Hypersensitivities

Article Number: 52 | Rating: Unrated | Last Updated: Tue, Jul 19, 2016 at 8:22 AM

Overview


The form allows entering new allergy/hypersensitivity information and creating reports and reviews of all collected data.

Applicable Staff


The ID Nurse can enter allergies and hypersensitivities. All users have viewing access.

Usage


The form is client specific (non-episodic) and allows adding the information to active and inactive episodes or both on the first screen. Open the form from a hyperlink in the sidebar of one of the other forms that is connected to it, from within a client record, or the menu path Avatar CWS/Other Chart Entry/Allergies and Hypersensitivities. The form opens as shown below:



When the **Update** button is clicked, the following multiple iteration table opens to enter new data or review/delete previous entries. For detailed information regarding using multiple iteration tables, see this article. Using a Form



After clicking **Save** in the table above, the view returns to the Allergies and Hypersensitivities screen of the form. Click the blue **Submit** button to save the info into the record. When allergy/hypersensitivity data is entered, it will be shown in the **Client Data Bar** on all screens.



Tips & FAQ

Any useful tips, and questions provided by users after roll out.

Technical Information

Menu Location (Avatar CWS/Other Chart Entry/Allergies and Hypersensitivities)

SQL Table(s)

RAD Envelope (if applicable)

Related Forms and Reports (Individual Information Sheet, Medical History Summary)

Posted by: [John Cookson](#) - Tue, Apr 26, 2016 at 6:47 PM. This article has been viewed 3216 times.

Online URL: <https://www.mercykb.com/article.php?id=52>